



Magnetic Stirrers “Agimatic” and “Agiman” analogue control

FOR SPEEDS FROM 60 TO 1600 rpm..
MAXIMUM STIR VOLUME: 10 LITRES.
PLATE MADE FROM AISI 304 STAINLESS STEEL WITH REMOVABLE SAFETY RING.
PLATE Ø 14.5 cm.



“Agimatic-S” without heating

CONTROL PANEL

1. Mains switch illuminated.
2. Electronic stir speed control.

Part No. **7000242**

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.



“Agimatic-N” with heating

ELECTRONIC DIGITAL CONTROL OF SPEED AND POWER TO THE HEATER.

CONTROL PANEL

1. Switch on button with indicator lamp for heating and for stirring speed.
2. Electronic stir speed control.
3. Electronic heater energy control from 0-100%, synchronised with an indicator lamp that functions when the plate is hot.
4. Illuminated indicator when the heater is “on”.

Part No. **7000243**

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.

Connection at the back to connect:

Electronic thermometer “Sensoterm II”



Electronic temperature regulator “Electemp”



“Agimatic-E” with heating

ELECTRONIC TEMPERATURE CONTROL OF SPEED AND POWER TO THE HEATER.
PRECISION ±2 °C.

CONTROL PANEL

1. Switch on button with indicator lamp for heating and for stirring speed.
2. Electronic stir speed control.
3. Electronic temperature control of the heater plate from 50 °C to 350 °C synchronised with an indicator lamp that functions when the plate is hot.
4. Illuminated heater indicator when the heater is on.

Part No. **7002431**


Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.

Connection at the back to connect:

Electronic thermometer “Sensoterm II”



Electronic temperature regulator “Electemp”



Magnetic stirrer “Agiman”

ELECTRONIC HOTPLATE DIGITAL CONTROL.
WITH A MECHANISM FOR ATTACHING HEATING MANTLES.

CONTROL PANEL

- Switch on button with indicator lamp for heating and for stirring speed.
Electronic stir speed control.
Electronic heater energy control from 0-100%, synchronised with an indicator lamp that goes on when the mantle is hot.
Illuminated heater indicator when the mantle is “on”.

Part No. **7000446**

Supplied with a PTFE coated stir bar of Ø 8 x 40 mm long.

Connection at the back to connect:

Electronic thermometer “Sensoterm II”



Electronic temperature regulator “Electemp”



ACCESSORIES

Attachable heating mantle for the “Agiman” magnetic stirrer.

Part No.	Capacity flask ml
3000447	250
3000448	500
3000449	1000
3000450	2000

Complements for Agimatic-S-N-E and Agiman.
See accessories on page 34 and 35.

