

# MAGNETIC STIRRERS



*A new generation of stirrer for the most exacting work*



## Magnetic Stirrer "Agimin" without heating

MICROPROCESSOR CONTROLLED, ELECTRONIC CONTROL, PORTABLE.

### FEATURES

Maximum stir volume 1 litre.  
Speed controllable from 40 to 2000 r.p.m. with storage of the ultimate speed setting.  
Double power source: 4 x 1.5V AA batteries or with a mains adapter.  
Moulded ABS plastic case.  
Excellent resistance to chemical agents.

### CONTROL PANEL

Luminous On/Off push button.  
Luminous r.p.m. push button.  
Increase/ decrease speed push button.  
Low battery indicator lamp.



### MODEL

Part No.	Height cm	Ø cm	Depth cm	Power W	Weight Kg
<b>7001731</b>	5	10	15	0.9	0.35

Comes complete with support bar and clamp for 3 electrodes (1 x 13 mm Ø and 2 x 16 mm Ø), 4 batteries 1.5V AA, mains adapters and 2 PTFE coated stir bars.



## Magnetic stirrer "Agimicro" without heating

### FEATURES

Maximum stir volume: 2 Litres.  
Adjustable stir speed from 40 to 1400 r.p.m.  
Durable moulded ABS plastic body.  
Excellent chemical resistance.

### CONTROL PANEL

Mains switch with luminous "on" indicator.  
Analogue control of speed in r.p.m.

### MODEL

Part No.	Height cm	Ø Plate cm	Power W	Weight Kg
<b>7001638</b>	5	12	2	0.5

Comes complete with a 8 Ø x 32 mm P.T.F.E. coated stir bar.

### ACCESSORIES

**Support bar and electrode clamp** (1 of 13 and 2 of 16 mm Ø). Used with pH and conductivity meters.  
Part No. **7001639**



## Powerful magnetic stirrer "Agimatic-H" without heating

FOR LOW VISCOSITY VOLUMES FROM 1 TO 25 LITRES AND 60 LITRES.  
USING REDUCED VOLUMES THE UNIT CAN STIR MEDIUM TO HIGH VISCOSITY LIQUIDS.

### Maximum stirring power

### FEATURES

Electronic speed control from 50 to 980 r.p.m. will maintain a constant stir speed irrespective of the liquid viscosity.  
Digital speed display.  
AISI 304 stainless steel top plate with an epoxy coated metal case.

### CONTROL PANEL

Main switch with a luminous "on" indicator.  
Electronic speed control.  
Digital speed display in r.p.m.

### MODEL

	Part No.	Speed control r.p.m.	Maximum stir volume in litres	Height/Width/Depth (exterior) cm	Power W	Weight Kg
<b>AGIMATIC-HS</b>	<b>7001005</b>	50 to 980	up to 25	12 28 34	35	8
<b>AGIMATIC-HL</b>	<b>7001006</b>	50 to 980	up to 60	14 44 38	40	12

Model HS comes complete with a 16 Ø x 127 mm stir bar and model HL with a 19 Ø x 150 mm stir bar.

