



Polyethylene beakers

Volumes available range from 60 ml to 1,000 ml.

All models come without lid, except the 1,000 ml beaker, which is available with or without lid.

For use with liquid and solid samples.

Manufactured in a flexible material for easy pouring.

code	volume ml	body Ø mm	mouth Ø mm	height mm	weight body g	case quantity	case weight	case volume
242801	60	40	50	48	5.0	1,100	5.95	0.034
242802	100	40	55	58	6.5	1,000	7.20	0.034
242803	230	55	70	78	15.0	700	10.82	0.042
242804	275	53	72	98	22.7	500	11.66	0.042
242805	500	70	90	108	38.7	200	8.50	0.045
242806	1,000	90	120	132	63.2	120	8.22	0.042
242811*	1,000	90	120	132	65.0	65	5.93	0.042

* With lid. Beaker and lid are supplied in the same bag.

Disposable cups

Disposable sample cups, for general sample collection.

Cups are supplied stacked.

code	volume ml	material	case quantity	case weight	case volume
202248	100	PP	50 x 100	7.40	0.080
202317	150	PS	30 x 100	8.32	0.080

