



Blue screw cap containers

CE (IVD)

With the same characteristics as codes **409702** (p. 137) and **409502** (p. 135) but in blue colour (containers and cap)

Ideal for photosensitive samples or when it is required not to show the content.

Ask for other labelling or for not labelled containers to our commercial department.

Minimum order quantity: 30.000 units.

code	description	case quantity	case weight	case volume	cases pallet
409511	60 ml container white label	600	8.00	0.082	24
409711	150 ml container white label	475	8.20	0.140	16




Graduated screw cap containers 57 x 116

CE (IVD)

Containers made of ultra clear polypropylene and caps made of polyethylene. Those containers have double moulded-in graduations, (in 20 ml increments up to 200 ml and in 1 oz increments up to 7 oz), and feature a frosted panel on their side for easy writing. Standard colour of caps: red.

code	max. capacity	case quantity	case weight	case volume	cases pallet
409902	240 ml	275	7.10	0.140	16



Dimensions:

A 56.5 mm; B 62.2 mm; C 115.5 mm;
Maximum volume 240 ml;
Recommended volume 210 ml.




Screw cap container 58 x 95 mm

CE (IVD)

Manufactured in polystyrene, cap made of white polyethylene.
Total volume: 220 ml

code	max. capacity	case quantity	case weight	case volume
309903	220 ml	240	9.6	0.108



Dimensions:

A 58 mm; B 62 mm; C 95 mm; D 22.