

Complete Mini Kits

Description	Titration method	Detection range	Additional reagent (number of pills)	Pack (u.)	Art. No.
Acid arsenic in drinking water	Colour	0,005-0,5mg/l As	No	1	644-400700
AF422 cyanuric acid	Turbidimetric	20-200mg/l cyanuric acid	No	1	644-414220
AF415 caustic alkalinity	Pills	20-500mg/l CaCO ₃	Art. No. 644-515501 (250 p)	1	644-414150
M alkalinity, AF 444	Sample drip	20-800mg/l CaCO ₃	No	1	644-414440
P alkalinity, AF 414	Pills	20-500mg/l CaCO ₃	Art. No. 644-515501 (250 p)	1	644-414140
AF413 total alkalinity	Pills	10-500mg/l CaCO ₃	Art. No. 644-515321 (250 p)	1	644-414130
AF417 quaternary ammonium	Pills	0-500mg/l AC active	No	1	644-414170
AF418 Chlorides	Pills	5-5.000mg/l CL	Art. No. 644-515131 (250 p)	1	644-414180
AF416 calcium Hardness	Pills	10-500mg/l CaCO ₃	Art. No. 644-515191 (250 p)	1	644-414160
AF446 calcium Hardness	Sample drip	20-800mg/l CaCO ₃	No	1	644-414460
AF423 hardness yes/no	Yes/No	4 or 8 or 20mg/l CaCO ₃	Art. No. 644-515361 (250 p)	1	644-414230
Total hardness, L AF445	Sample drip	20-800mg/l CaCO ₃	No	1	644-414450
Total hardness, P AF424	Pills	5-500mg/l CaCO ₃	Art. No. 644-515161 (250 p)	1	644-414240
Hardness, low range, AF425	Pills	1-50mg/l CaCO ₃	Art. No. 644-515171 (250 p)	1	644-414250
Hardness, very low range AF426	Pills	1-10mg/l CaCO ₃	Art. No. 644-515351 (250 p)	1	644-414260
AF411 organic phosphonate	Pills	1-20mg/l O-P	Art. No. 644-465351 and 644-512961 (250 p)	1	644-414110
AF427 Nitrites	Pills	70-1.500mg/l NaNO ₂	Art. No. 644-515201 and 644-515211 (250 p)	1	644-414270
AF410 acid cleaner	Pills	0,75 a 10%	No	1	644-414100
Sulfates, low level, AF432	Pills	20-200mg/l Na ₂ SO ₄	Art. No. 644-515221 and 644-515231 (250 p)	1	644-414320
Sulfates, AF431	Turbidimetric	40-4.000mg/l SO ₄	Art. No. 644-515451 (250 p)	1	644-414310
Sulfite, high level, AF435	Pills	20-500mg/l Na ₂ SO ₃	Art. No. 644-515271 and 644-515281 (250 p)	1	644-414350
Sulfite, low level, AF434	Pills	2-50mg/l Na ₂ SO ₃	Art. No. 644-515271 and 644-515331 (250 p)	1	644-414340
Tannins, AF436 index	Pills	2-20mg/l units	Art. No. 644-503511 (250 p)	1	644-414360

WATER ANALYSIS: Test strips



Semi-quantitative test strips SCHARLAU

Scharlau's semi-quantitative test strips enable an easy and quick determination of different parameters in a solution.

Description	Measuring values (ppm)	Pack (u.)	Art. No.
Nitrite indicator strips. 100 strips per package	Readings at 0-0,5-1-5-10-25 ppm	1	TP002500NI
Nitrate indicator strips. 100 strips per package	Readings at 0-10-25-50-100-250-500 ppm	1	TP050000NA
Peracetic acid indicator strips. 100 strips per package	Readings at 0-5-10-20-30-50 ppm	1	TP005000PA
Peracetic acid indicator strips. 100 strips per package	Readings at 0-100-250-500 ppm	1	TP050000PA
Peroxide indicator strips. 100 strips per vial	Readings at 0-1-3-10-30-100 ppm	1	TP010000PX
Peroxide indicator strips. 100 strips per vial	Readings at 0-100-200-400 ppm	1	TP040000PX

For further information, please contact scharlab@scharlab.com