



Blue Nitrile Disposable Xtra lite powder-free gloves SEMPERGUARD®

Resistant to chemicals, acids and non ketone solvents.
Larger pack, less packaging waste.
Ambidextrous, with rolled cuffs, impermeable.
Free from allergenic latex proteins
0,07mm wall thickness. Category III. Conformity to the norms EN 374 and EN 455. Suitable for food contact and AQL 1,5.

Size	Length (mm)	Pack (u.)	Art. No.
S	240	200	266SP0060S
M	240	200	266SP0060M
L	240	200	266SP0060L
XL	240	180	266SP006XL



Confort Blue Nitrile Disposable powder-free gloves SEMPERGUARD®

CE category III protective glove for complex risks
Soft glove with excellent chemical resistance to skin irritants.
High elasticity similar to latex products to guarantee comfort.
Protection against chemical splashing including the processing of fatty foods.
Free from phthalates/softeners and allergenic latex proteins.
Category III. Conformity to the norms EN 374 and EN 455. Suitable for food contact and AQL 1,5.

Size	Length (mm)	Pack (u.)	Art. No.
S	240	100	266SP6080S
M	240	100	266SP6080M
L	240	100	266SP6080L
XL	240	90	266SP608XL



XPert Blue Nitrile Disposable powder-free gloves SEMPERGUARD®

Resistant to chemicals, acids and non ketone solvents.
Ambidextrous. Impermeable. Powder-free to avoid allergies.
Category III. Conformity to the norms EN 374. Suitable for food contact and AQL 1,5.

Size	Pack (u.)	Art. No.
S	100	2666080SPS
M	100	2666080SPM
L	100	2666080SPL
XL	90	2666080SXL

Visit our new web www.scharlab.com

PERSONAL PROTECTIVE EQUIPMENT