



SURFACE CONTROL

2•4 Incubators - ATL® mini EM4

BENEFITS

- ❑ Economically priced
- ❑ Convenient and rational:
 - Removable lid with handle, for fast top loading and removal of slides and plates.
 - Easy to clean.
 - Rack for 23 contact slides and thermometer are included as a standard.

FIELDS OF APPLICATION

- ❑ Production, on-site analyses:
 - for the incubation of **contact slides or plates** after surface sampling.
 - for the incubation of contact plates after environmental sampling of air.
- ❑ In the microbiology laboratory:
 - for incubation of small amounts of Petri dishes.

TECHNICAL CHARACTERISTICS

- ❑ Adjusting screw for temperature setting (from room temperature +5° C up to 44°C).
- ❑ Volume: 4 litres.
- ❑ External dimensions: 23 x 23 x h. 18 cm.
- ❑ Internal dimensions: 18 x 18 x h. 12.5 cm.
- ❑ Weight <1.5 kg.

For Contact Slides

607.028053.01 ATL® Mini EM 4 Incubator; 4 litres; **includes rack** for 23 slides and thermometer



For Contact Plates

607.028053.11 ATL® Mini EM 4 Incubator; 4 litres; **includes bottom shelf** and thermometer



Accessory

607.028053.10 Shelf

Slide Racks

ATL® Transport Rack, stainless steel, with folding handle:

- Holds 24 slides.
- Electropolished.
- Size 245 x 155 x h. 85 mm.
- Effective grid size: 30 mm.

157.048500.02 ATL® Transport Rack for **24 contact slides**



ATL® Rack, for use with MINI EM4 Incubator

607.028053.02 Rack with centred handle for 24 slides

