



EUROTUBO® 50 ml conical tubes

CE (IVD)

Tubes made of transparent polypropylene, conical bottom, suitable for centrifugation tests in immunology, microbiology, etc. Perfect to analyse Koch bacillus. Continuous thread, external moulded graduations in 5 ml increments. Blue cap in high density polyethylene with hermetical closure thanks to its internal elastic obturation. Skirted and non skirted versions, sterilised by radiation and non sterilised versions.

Sterile codes are either supplied in individual bags (codes **429926**, **429927**), in bags of 100 units (codes **429926.10**, **429927.10**), or in bags of 25 units (codes **429926.25**, **429927.25**) and resist up to **7,000 xg**. Autoclave not recommended as they have been ionised. Non sterile codes withstand up to **12,000 xg** and are autoclavable (**121 °C**) with the cap loose on the thread and not tightly fitted.

Every version is supplied capped, excepting code **429900SP**, which is supplied uncapped, with the cap in a separated bag.

code		description	external cap Ø mm	internal tube Ø mm	external tube Ø mm	length with cap mm	case quantity	case weight	case volume	cases per pallet
429900		non skirted	34.4	27.2	29.5	117.5	500	7.9	0.075	20
429900SP		non skirted, unscrewed cap	34.4	27.2	29.5	117.5	500	7.9	0.082	20
429901		skirted	34.4	27.2	29.5	117.5	500	8.5	0.082	20
429926	STERILE	non skirted, individually wrapped	34.4	27.2	29.5	117.5	500	8.2	0.082	20
429926.25	STERILE	non skirted, bag 25 units	34.4	27.2	29.5	117.5	20x25	8.2	0.082	20
429926.10	STERILE	non skirted, bag 100 units	34.4	27.2	29.5	117.5	5x100	7.9	0.082	20
429927	STERILE	skirted, individually wrapped	34.4	27.2	29.5	117.5	500	8.8	0.082	20
429927.25	STERILE	skirted, bag 25 units	34.4	27.2	29.5	117.5	20x25	8.8	0.082	20
429927.10	STERILE	skirted, bag 100 units	34.4	27.2	29.5	117.5	5x100	8.5	0.082	20



Please find racks for these tubes in Chapter 9. **Sample Storage**



POLYPROPYLENE RACKS:

- W-018
- 19568
- 19570

METAL RACKS:

- R-292 • R-293
- R-281 • R-282
- R-283

