



Large filter

Description	Pack (u.)	Art. No.
Small filter for containers <5 liters. 3 months shelf life	1	1368107911
Small filters for containers <5 liters. 3 months shelf life	4	1368190911
Medium filter for containers <20 liters. 6 months shelf life	1	1368107914
Medium filters for containers <20 liters. 6 months shelf life	3	1368190914
Large filter for containers >20 liters. 9 months shelf life	1	1368107615

## WASTE: Waste containers



### Biosafety containers

Made of rigid polypropylene. Shock, puncture and solvent resistant. Incinerable, without residual waste. Suitable for solid waste from groups II and III. Manufactured according to ISO 9001 and approved according to ADR for transport.

Description	Dimensions (mm)	Capacity (l)	Pack (u.)	Art. No.
Oval aperture and for syringes	107x90x190	1	1	027-011000
Five types of openings for items of different sizes	195x120x170	2	1	027-012000
Five types of openings for items of different sizes	175x175x248	4	1	027-040004
Five types of openings for items of different sizes	175x175x382	7	1	027-040007



### Biosafety containers with handle

Made of virgin totally incinerable materials, without residual waste. Suitable for solid waste types II and III. Complies with almost all European and international regulations (ADR;ONU).

Description	Height (mm)	Capacity (l)	Pack (u.)	Art. No.
Not autoclavable	162	2	1	027-022000
Autoclavable	172	4	1	027-024000
Autoclavable	201	6	1	027-026000
With full or partial opening of the lid	255	11	1	027-030011

For further information, please contact [scharlab@scharlab.com](mailto:scharlab@scharlab.com)