


### CIRCULAR SUPPORT PLATFORMS WITH BOSSHEAD

To support flat bottom vessels and instruments.  
Chrome plated bosshead.

Part No.	Plate Ø
6000879	80 mm.
6000880	150 mm.



### OPEN RING BOSSHEAD SUPPORTS

Made of galvanised steel.


Part No.	Ring exterior Ø	Usable length A
7000506	60 mm.	100 mm.
7000508	80 mm.	130 mm.
7005010	100 mm.	160 mm.
7005012	120 mm.	170 mm.
7005014	140 mm.	180 mm.



### CLOSED RING BOSSHEAD SUPPORTS

Made of galvanised steel.


Part No.	Ring exterior Ø	Usable length A
7000536	60 mm.	100 mm.
7000538	80 mm.	130 mm.
7005310	100 mm.	160 mm.
7005312	120 mm.	170 mm.
7005314	140 mm.	180 mm.



### OPEN RING RETORT SUPPORTS

Made of galvanised steel.


Part No.	Ring exterior Ø	Usable length A
1000206	60 mm.	180 mm.
1000208	80 mm.	195 mm.
1002010	100 mm.	200 mm.
1002012	120 mm.	215 mm.
1002014	140 mm.	225 mm.



### CLOSED RING RETORT SUPPORTS


Made from galvanised steel.

Part No.	Ring exterior Ø	Usable length A
1000236	60 mm.	180 mm.
1000238	80 mm.	195 mm.
1002310	100 mm.	200 mm.
1002312	120 mm.	215 mm.
1002314	140 mm.	225 mm.



### GAUZE IRON WIRE WITH CERAMIC FIBRE

Part No.	Side
1087102	100 mm Side.
1087122	120 mm Side.
1087142	140 mm Side.
1087162	160 mm Side.
1087182	180 mm Side.
1087202	200 mm Side.




### CERAMIC PLATES "NEOCERAM"

Resistant to temperatures from -200 °C to 700 °C.

Part No.	Plate Size
1000684	120 x 120 x 5 mm.
1000685	140 x 140 x 5 mm.
1000686	160 x 160 x 5 mm.

### Heat conductivity



Resistant to temperatures from -200 °C to 700 °C.

Excellent heat conductor.

20% better heat transfer compared to a conventional gauze ceramic.

High chemical resistance to aggressive compounds.

Easy to clean and long lasting.



### TRIPOD STAND FOR CERAMIC PLATES

With adjustable sides to prevent the plate from moving or tilting.

Total height 210 mm.

Part No.	For plates of
1000687	120 x 120 mm.
1000688	140 x 140 mm.
1000689	160 x 160 mm.