



Plastic bucket

Plastic bucket of high density polyethylene, with lid.

Capacity (Kg)	Pack (u.)	Art. No.
5	1	0559PI5000
10	1	05591000CP
25	1	05592500CP



Flasks

Cylindrical flasks, made of high density polyethylene with knurled screw cap and stopper.

Ø int.xH (mm)	Cap. (ml)	Pack (u.)	Art. No.
50x94	250	50	425-001565
70x107	500	50	425-001566
85x128	1000	20	425-001567
54x100	250	125	027-202802
67x131	500	120	027-202803
86x157	1000	50	027-202818
60x124	400	130	027-202814
86x103	500	100	027-202821
86x173	1.250	50	027-202804

Sterile single bag

Ø int.xH (mm)	Cap. (ml)	Pack (u.)	Art. No.
54x100	250	125	00202802.1
67x131	500	120	00202823.1
86x157	1000	50	00202828.G

Single bag

Ø int.xH (mm)	Capacity (ml)	Pack (u.)	Art. No.
54x100	250	125	027202802B
67x131	500	120	027-202823
86x157	1000	50	027-202828



Flasks GOSSELIN SAS

Cylindrical flasks (Dutchess type) of high density PE with red caps.
Different closure systems available.
All caps are knurled for an easy opening even with gloves.

Ø int.xH (mm)	Cap type	Cap. (ml)	Pack (u.)	Art. No.
60x25	Cap with wad	250	145	630LR25004
80x136	Cap with wad	500	140	630LR50006
101x174	Cap with wad	1000	68	00LR100004
60x25	Cap with seal	250	145	630LR25001
80x136	Cap with seal	500	140	630LR50001
101x174	Cap with seal	1000	68	00LR100001
80x136	Cap with wad. Sterile	500	140	630LR50004
101x174	Cap with wad. Sterile	1000	68	00LR100008
60x25	Cap with seal. Sterile	250	145	630LR25006
80x136	Cap with seal. Sterile	500	140	630LR50003
101x174	Cap with seal. Sterile	1000	68	00LR100002

For further information, please contact scharlab@scharlab.com