



Graduated wide neck bottles

KARTELL

Graduated bottle for liquids and powders. Shoulder and cap with holes for safety seal. Permanent graduation printed in relief externally. Will withstand considerable pressure before rupture.

Ø int.xH (mm)	Cap. (ml)	Pack (u.)	Art. No.
61x126	250	50	425-000410
75x158	500	20	425-000411
93x200	1000	20	425-000412



Square bottles

DELALAB

Square flasks with star-shaped cap and high-density polyethylene seal (PE). The square bottle is space-saving. Suitable for long term sample storage. Graduation of the flasks permits filling without the need for pipettes or other utensils

Cap. (ml)	Pack (u.)	Art. No.
100	600	027-444610
250	300	027-444611
500	150	027-444612
1000	80	027-444613



Square bottles

GOSSELIN SAS

Square flasks of PEHD, graduated.

Non-sterile flasks

Dim. WxDxH (mm)	Cap type	Diameter (mm)	Cap. (ml)	Pack (u.)	Art. No.
57x57x109	Cap with wad	37	250	210	630LC25010
76x76x122	Cap with wad	58	500	175	630LC50007
95x76x185	Cap with wad	58	1000	90	00LC100007
109x96x230	Cap with wad	58	2000	50	00LC200005
57x57x109	Cap with seal	37	250	210	630LC25001
76x76x122	Cap with seal	58	500	175	630LC50001
95x76x185	Cap with seal	58	1000	90	00LC100001
109x96x230	Cap with seal	58	2000	50	00LC200001

Sterile flasks

Dim. WxDxH (mm)	Cap type	Diameter (mm)	Cap. (ml)	Pack (u.)	Art. No.
57x57x109	Cap with wad	37	250	210	630LC25006
76x76x122	Cap with wad	58	500	175	630LC50011
95x76x185	Cap with wad	58	1000	90	00LC100013
57x57x109	Cap with seal	37	250	210	630LC25002
76x76x122	Cap with seal	58	500	175	630LC50002
95x76x185	Cap with seal	58	1000	90	00LC100002
109x96x230	Cap with seal	58	2000	50	00LC200003

Visit our new web www.scharlab.com