

Zirconium lids for crucible

J.P. SELECTA

Ø (mm)	Pack (u.)	Referencia
33	1	0196000567
47	1	0196000568
52	1	0196000569
59	1	0196000570

CRUCIBLES: Porcelain crucibles

Alsint 99,7 ceramics crucibles

HALDENWANGER

Temperature resistance up to 1.700°C.

Combustion boats



Combustion boat

Length (mm)	Width (mm)	Height (mm)	Pack (u.)	Art. No.
50	20	20	1	011-CAJA50
50	25	20	1	011-CAJA51
100	32	28	1	011-CAJ100
150	65	19	1	011-CAJ150

Conical, squat shape crucibles



Wide form crucible

Ø sup. / Ø inf. (mm)	Height (mm)	Cap. (ml)	Pack (u.)	Art. No.
41/18	37	25	1	031-00000B
48/20	41	40	1	011-00001B
54/24	50	60	1	011-00002B
60/26	50	80	1	011-00003B

Conical, tall shape crucibles



Ø sup. / Ø inf. (mm)	Height (mm)	Cap. (ml)	Pack (u.)	Art. No.
33/18	40	20	1	011-00001C
38/21	47	30	1	011-00002C
42/25	54	45	1	011-00003C
50/27	65	80	1	011-00004C

Visit our new web www.scharlab.com