

Typ 2

### Messbecher graduiert, konische Form

Graduierung 100 ml

### Measuring jug graduated, conical shape

Graduation 100 ml


Material	C	D	MI	Typ	Artikel Nr.
18/10-Stahl	120	125	1000	1	8640
18/10-Stahl	145	125	2000	1	8642
18/10-Stahl	125	90	500	1	8650
18/10-Stahl	140	125	1000	1	8651
18/10-Stahl	170	140	2000	1	8652
18/10-Stahl	115	115	500	2	8660
18/10-Stahl	160	140	1000	2	8661
18/10-Stahl	160	150	1500	2	8662
18/10-Stahl	185	160	2000	2	8663

### Pipettenbüchse, rund

Silikon im Boden u. Deckel, bis 205 °C sterilisierbar

### Pipette box, round


Silicone covered bottom and lid, sterilizable up to 205 °C



Material	C	c	D	Artikel Nr.
18/10-Stahl	390	382	50	8980
18/10-Stahl	425	417	50	8981
18/10-Stahl	380	372	65	8982
18/10-Stahl	450	442	65	8983
18/10-Stahl	380	372	80	8984
18/10-Stahl	500	492	80	8985



Lieferung ohne Pipetten / delivery without pipettes



### Pipettenbüchse, eckig

Silikon im Boden und Deckel, bis 205 °C sterilisierbar

### Pipette box, square

Silicone covered bottom and lid, sterilizable up to 205 °C

Material	C	c	D	Artikel Nr.
18/10-Stahl	210	202	70 x 70	8870
18/10-Stahl	290	282	70 x 70	8871
18/10-Stahl	340	332	70 x 70	8872
18/10-Stahl	435	427	70 x 70	8873
18/10-Stahl	490	482	70 x 70	8874

### Sterilisierbehälter für Petrischalen

Ø 100 mm (8541)

Für 10 Schalen

### Sterilizing box for Petri-dishes

Ø 100 mm (8541)

For 10 dishes

Material	C	D	Typ	Artikel Nr.
18/10-Stahl	250	120	10 x	8961



Petrischalen siehe S. 60 /  
petri dishes see p. 60

