




Open ring bosshead supports made of galvanized steel

Ø ext. ring (mm)	Usable length (mm)	Pack (u.)	Art. No.
60	100	1	0197000506
80	130	1	0197000508
100	160	1	0197005010
120	170	1	0197005012
140	180	1	0197005014



Open ring retort supports made of galvanized steel

Ø ext. ring (mm)	Usable length (mm)	Pack (u.)	Art. No.
60	180	1	0191000206
80	195	1	0191000208
100	200	1	0191002010
120	215	1	0191002012
140	225	1	0191002014



Support rings

Stainless steel 18/10. With and without bosshead.

Ø int. ring (mm)	Bosshead	Pack (u.)	Art. No.
70	With clamp	1	5501000BOC
100	With clamp	1	5502000BOC
120	With clamp	1	5503000BOC
140	With clamp	1	5504000BOC
70	Without clamp	1	5501M00BOC
100	Without clamp	1	5502M00BOC
120	Without clamp	1	5503M00BOC
140	Without clamp	1	5504M00BOC

INSTRUMENTAL: Platforms



Aluminium anodised Lab-Jacks

Aluminium anodised. DIN 12897. Blue. Smaller models up to 240x240mm have a simple adjusting wheel, which can be operated by hand. With the bigger models from 300x300mm on we recommend the use of a ratchet, which connects with a special adjusting wheel.

NOTE: Please order ratchet separately.

Description	Size (mm)	Elevator height (mm)	Dynamic max. capacity (Kg)	Static max. capacity (Kg)	Pack (u.)	Art. No.
Lab-Jack	100x100	55 to 120	5	10	1	1101600BOC
Lab-Jack	160x130	60 to 275	7	30	1	1102100BOC
Lab-Jack	200x200	60 to 275	7	30	1	1103100BOC
Lab-Jack	240x240	60 to 275	7	30	1	1104100BOC
Lab-Jack	300x300	130 to 470	12	60	1	1108100BOC
Ratchet	190	-	-	-	1	1108800BOC