




Set for students, 7 pieces

Simple type:
 Scalpel handle No. 3.
 5 pieces each of non-sterile scalpel blades Fig. 10 and 15.
 Preparations scissors.
 Fine, sharp forceps.
 Lancet dissection needle.
 Dissection needle.

Type	Pack (u.)	Art. No.
Set for students, 7 pieces	1	HSO-130-00

INSTRUMENTAL: Spatulas



Chemical spoon, titanium, 180mm long

Width (mm)	Spoon (mm)	Pack (u.)	Art. No.
22	28x40	1	3920000BOC



Double ended spoons, steel 18 / 10

Length (mm)	Spatula (mm)	Pack (u.)	Art. No.
120	22x30 and 17x23	1	3460000BOC
150	22x30 and 17x23	1	3461000BOC
180	29x40 and 22x30	1	3462000BOC
210	29x40 and 22x30	1	0005150021
250	29x48 and 29x40	1	3464000BOC
300	35x48 and 29x40	1	3465000BOC



Double spatula, titanium

Size (mm)	Ø (mm)	Pack (u.)	Art. No.
11x210	4	1	3910000BOC

Visit our new web www.scharlab.com