

**MICROBIOLOGY: Petri dishes and accessories**
**Petri dishes**  
**GOSSELIN SAS**

Gosselin Petri dishes are manufactured following a totally automated process with strict quality controls. Made of high transparency PS, using a lateral injection process which confers the ideal consistency and transparency for microbiological applications. They do not deform at 60°C. Excellent mechanical resistance. Excellent optical qualities. Available with or without vents, making air circulation easier and avoiding condensation.

Different sterility levels:

- Aseptic: manufactured under strict aseptic conditions, using a "Clean Room" technology, production and packaging line process. The plastic injection temperature is above 250°C. This quality is recommended for use in bacteriological applications where contamination risk is almost nonexistent.
- Sterile: dishes undergo an irradiation dose of 10Kgray. Each box comes with a certificate accrediting said irradiation.



63090143NA



630-902CNA



630-55143A



630-65143A



630-14203A

Description	Ø (mm)	Height (mm)	Amount/Bag	Pack (u.)	Art. No.
Aseptic ISO 5, 3 vents	90	14,2	Bags 33	825	63090143NA
Sterile, 3 vents	90	14,2	Bags 33	825	63090143NE
Aseptic ISO 5, without vents	90	14,2	Bags 33	825	63090140NA
Sterile, without vents	90	14,2	Bags 33	825	63090140NE
Aseptic ISO 5, 3 vents	90	16,2	Bags 25	700	630-90163A
Sterile, 3 vents	90	16,2	Bags 25	700	630-90163E
Aseptic ISO 5, 3 vents, 2 compartments	90	14,2	Bags 33	825	630-902CNA
Sterile, 3 vents, 2 compartments	90	14,2	Bags 33	825	630-902CNE
Aseptic ISO 5, 6 vents	55	14,2	Bags 15	1620	630-55143A
Sterile, 6 vents	55	14,2	Bags 15	1620	630-55143E
Aseptic ISO 5, 3 vents	65 (Contact)	14,5	Bags 20	720	630-65163A
Sterile, 3 vents	65 (Contact)	14,5	Bags 20	720	630-65163E
Aseptic ISO 5, 4 vents with safety lock "Grip-Lid" and convex base	65 (Contact)	14,5	Bags 20	720	630-65143A
Sterile, 4 vents with safety lock "Grip-Lid" and convex base	65 (Contact)	14,5	Bags 20	720	630-65143E
Aseptic ISO 5, 3 vents	140	20,4	Bags 11	176	630-14203A
Sterile, 3 vents	140	20,4	Bags 11	176	630-14203E
Aseptic ISO 5, 4 vents, square	120	15,8	Bags 14	252	630-12154A
Sterile, 4 vents, square	120	15,8	Bags 14	252	630-12154E


**Stand for dishes**

Made of polypropylene, it increases the incubator's capacity thanks to its stackable and secure design.

Designed to contain five 90mm petri dishes, or fifteen 60mm dishes. Provided with two areas for labeling. Non autoclavable.

Description	Dim. WxHxD (cm)	Pack (u.)	Art. No.
PP stand for incubation	25,1x23,7x3,5	3	027-0H-615

Visit our new web [www.scharlab.com](http://www.scharlab.com)