


**Pipetus® automatic pipette filler**  
**HIRSCHMANN**

Rechargeable battery operated pipetting filler. All available pipettes are suitable from 0,1 to 200ml. Speed switch selection between 5 suction and dispensing speed possibilities. More than 50h work without charging. Without noise and vibrations. Very light. With safety valve which prevents liquid from entering the unit. Finest meniscus adjustment by superior valve technology. Hydrophobic PTFE 0,2µm filters. Autoclavable pipetting filler except from PTFE filters. Silicone adapter ensures the seating of all pipettes. Provided complete with: 2 filters, safety valve and charger. With integrated colour display which shows status of battery and buzzer when battery reaches under 10%.

Description	Pack (u.)	Art. No.
Pipetting filler Pipetus® (Complete)	1	6709907200


**Pipetus®-Standard automatic pipette filler**  
**HIRSCHMANN**

For working explosive and/or sterile environments. Motoric pipette filler for volumes between 0,1ml-200ml. Powered by an external electric pump which makes it ideal for working in potentially explosive environments, with toxic reagents or in clean room technology. Ideal for applications where sterility is an essential requirement, due to its pump can incorporate an extremely effective filter, which can be integrated on both air directions. Complete with electric pump, 2 filters and an extra safety valve.

Description	Pack (u.)	Art. No.
Pipeteador Pipetus®-Standard with pump (complete)	1	6709903000
Pipeteador Pipetus®-Standard with pump and external filter	1	6709904000

**Accessories for Pipetus® line**  
**HIRSCHMANN**

Description	Pack (u.)	Art. No.
Filter 0,2µm	5	6709903704
Holder for all Pipetus® models	1	6709905700


**Automatic pipettor Seropet**  
**KARTELL**

Battery life: 8h. Autoclavable pieces: filters, pipette holder and nosepiece. Filters: hydrophobic PTFE 0,2µm. Compatible pipettes: glass and plastic from 1 to 10ml. High aspiration speed: 25ml in 4 seconds. Great valve sensitivity which allows high working precision with low volume pipettes. LED adviser when low battery Includes: NiMH battery, charger, hydrophobic PTFE 0,2µm filter and adhesive holder.

Description	Pack (u.)	Art. No.
Seropet 1-100ml	1	425-129510
2 PTFE autoclavable filters 0,2µm and 0,45µm	1	425-001360

## Glass pipettes in PIPETTES

