

PIPETTES: Serological pipettes

**Serological pipettes
DELTALAB**

Single use, made of crystal polystyrene. Accuracy > 98%. Sterilised by radiation. Made in one piece for an optimum precision, except Art. No. 90043 (50ml). No liquid retention: as tip is not sealed, liquid is not held at the joint. Total dispensing guaranteed. Non pyrogenic (ISO 10993-11), non cytotoxic (ISO 10993-5) and non haemolytic (ISO 10993-4). DNase and RNase free. Guaranteed features for 5 years. Each pack is delivered with sterility and conformity certificate.

Capacity (ml)	Tip	Total cap. (ml)	Presentation	Pack (u.)	Art. No.
1	Tape end	1,3	Unit bag	500	027-900030
1	Tape end	1,3	Bags 25 units	1.000	027-900031
1	Open tip	1,3	Unit bag	500	027-900130
1	Open tip	1,3	Bags 25 units	1.000	027-900131
2	Tape end	2,4	Unit bag	500	027-900032
2	Tape end	2,4	Bags 25 units	1.000	027-900033
5	Tape end	8	Unit bag	250	027-900034
5	Tape end	8	Bags 25 units	375	027-900038
5	Wide tip	8	Unit bag	125	027900144N
5	Open tip	8	Unit bag	125	027900134N
10	Tape end	13	Unit bag	200	027-900036
10	Tape end	13	Bags 25 units	250	027-900037
10	Wide tip	13	Bags 25 units	250	027900146N
10	Open tip	13	Unit bag	125	027900136N
10	Open tip	13	Bags 25 units	250	027-900137
25	Tape end	35	Unit bag	100	027-900041
50	Tape end	55	Unit bag	200	027-900043

PIPETTES: Volumetric pipettes

**Volumetric pipettes, one mark, class AS, serialized and certified, soda-lime glass
SCHARLAU**

Capacity (ml)	Pack (u.)	Art. No.
1	5	073-991676
2	5	073-991677
5	5	073-991678
10	5	073-991679
20	5	073-991681
25	5	073-991682
50	2	073-991683


**Volumetric pipettes, two marks, class AS, serialized and certified, soda-lime glass
SCHARLAU**

Capacity (ml)	Pack (u.)	Art. No.
1	5	073-991696
2	5	073-991697
5	5	073-991698
10	5	073-991699
20	5	073-991701
25	5	073-991702
50	2	073-991703