

**Measuring scoops of polyethylene**
**VITLAB**

High density polyethylene (PEHD). Natural colouring, industrial version. Suitable for contact with foodstuff in accordance with Regulation (EC) No. 10 / 2011.



Length (mm)	Capacity (ml)	Pack (u.)	Art. No.
115	15	1	309-040093
135	25	1	309-040193
185	65	1	309-040293
215	110	1	309-040393
250	150	1	309-040493
310	350	1	309-040593
350	750	1	309-040693
400	1.250	1	309-040793

**Measuring scoops of polypropylene**
**VITLAB**

Polypropylene. White, flat bottoms and tops, volume-marked on handle. Suitable for contact with foodstuff in accordance with Regulation (EC) No. 10/2011.



Length (mm)	Capacity (ml)	Pack (u.)	Art. No.
60	2	1	0005201052
82	5	1	0005423770
100	10	1	0005201055
135	25	1	309-039494
160	50	1	309-039594
200	100	1	309-039694
260	250	1	309-039794
315	500	1	309-039894
385	1.000	1	309-039994

**Aluminium scoops**
**WILESCO**

Aluminium.



Length (mm)	Pack (u.)	Art. No.
135	1	4279135WIL
185	1	4279185WIL
210	1	4279210WIL
250	1	4279250WIL
265	1	4279265WIL
310	1	4279310WIL
350	1	4279350WIL
400	1	4279400WIL
430	1	4279430WIL

Visit our new web [www.scharlab.com](http://www.scharlab.com)