

Test tube racks


J.P. SELECTA

For tubes up to 13mm in diameter



For tube Ø up to (mm)	Number of tubes	Dim. WxDxH (mm)	Made of	Pack (u.)	Art. No.
13	12 (4x3)	65x50x65	Plasticized wire	1	0191001212
13	12 (4x3)	65x50x65	Stainless steel wire	1	0191025500

Up to 12 tubes of 16 or 25mm diameter



For tube Ø up to (mm)	Number of tubes	Dim. WxDxH (mm)	Made of	Pack (u.)	Art. No.
16	12 (6x2)	120x45x80	Plasticized wire	1	0191026112
25	12 (6x2)	180x65x110	Plasticized wire	1	0191026412
16	12 (6x2)	120x45x80	Stainless steel wire	1	0191025512
25	12 (6x2)	180x65x110	Stainless steel wire	1	0191025812

Up to 24 tubes of 16 or 25mm diameter



For tube Ø up to (mm)	Number of tubes	Dim. WxDxH (mm)	Made of	Pack (u.)	Art. No.
16	24 (8x3)	160x65x80	Plasticized wire	1	0191026124
25	24 (8x3)	235x95x110	Plasticized wire	1	0191026424
16	24 (8x3)	160x65x80	Stainless steel wire	1	0191025524
25	24 (8x3)	235x95x110	Stainless steel wire	1	0191025824

Up to 48 tubes of 16 or 25mm diameter



For tube Ø up to (mm)	Number of tubes	Dim. WxDxH (mm)	Made of	Pack (u.)	Art. No.
16	48 (8x6)	160x125x80	Plasticized wire	1	0191026148
25	48 (8x6)	235x180x110	Plasticized wire	1	0191026448
16	48 (8x6)	160x125x80	Stainless steel wire	1	0191025548
25	48 (8x6)	235x180x110	Stainless steel wire	1	0191025848

Up to 50 tubes of 13 or 18mm diameter



For tube Ø up to (mm)	Number of tubes	Dim. WxDxH (mm)	Made of	Pack (u.)	Art. No.
13	50 (10x5)	162x84x35	Stainless steel wire	1	0191000603
13	50 (10x5)	162x84x60	Stainless steel wire	1	0191000605
18	50 (10x5)	210x110x40	Stainless steel wire	1	0191000607
18	50 (10x5)	210x110x60	Stainless steel wire	1	0191000609
18	50 (10x5)	210x110x85	Stainless steel wire	1	0191000611