



Sterile culture tubes in polystyrene DELTA LAB

Tubes supplied with either a two position ribbed polyethylene cap, which can be left loose for aerobic work or sealed for anaerobic cultures. Can withstand up to 1.400g centrifugation and 70°C.

Dim. ØxH (mm)	Pack (u.)	Art. No.
12x75	1.000	027-300807
17x100	1.000	027-300808



Straigh container 60ml polystyrene and screw cap GOSSELIN SAS

Description	Dim. ØxH (mm)	Pack (u.)	Art. No.
PS	33x70	700	00TP35-002
PS. Sterile	33x70	700	00TP35-103



30ml universal polystyrene containers DELTA LAB

Description	Dim. ØxH (mm)	Pack (u.)	Art. No.
PS	30,8x92,5	400	027-309202
PS. Unit bag	30,8x92,5	400	027-309222
PS. Unit bag. Sterile	30,8x92,5	400	00309222.0



Durham tubes SCHARLAU

Description	Bottom	Cap. (ml)	Dim. ØxLength (mm)	Material	Tubes/box	Pack (u.)	Art. No.
Durham (bell) tube	Round	0,6	6x20	Soda	100	1	124-00620C
Durham (bell) tube	Round	0,9	6x30	Soda	100	1	124-00630C
Durham (bell) tube	Round	1,1	6x40	Soda	100	1	124-00640C
Durham (bell) tube	Round	1,1	7x30	Soda	100	1	124-00730C
Durham (bell) tube	Round	1,4	7x40	Soda	100	1	124-00740C
Durham (bell) tube	Round	1,5	8x30	Soda	100	1	124-00830C
Amber glass tube	Round	4	11x50	Amber soda	100	1	073-003654

For further information, please contact scharlab@scharlab.com