

TUBES: Brushes



Brushes

Brushes made of natural hair for glass material cleaning. Metallic and flexible handle.

For big flasks

Total length (mm)	Usable length (mm)	Ø (mm)	Pack (u.)	Art. No.
620	130	100	5	233-402-14
550	120	90	5	233-402-13
500	100	75	5	233-402-12
370	100	70	5	233-0402-7

For middle flasks

Total length (mm)	Usable length (mm)	Ø (mm)	Pack (u.)	Art. No.
320	90	60	5	233-0402-6
360	100	45	5	233-0402-4
250	75	50	5	233-0402-5

For small flasks

Total length (mm)	Usable length (mm)	Ø (mm)	Pack (u.)	Art. No.
270	95	35	5	233-0402-3

For burettes

Total length (mm)	Usable length (mm)	Ø (mm)	Pack (u.)	Art. No.
670	100	23	5	233-402-10

For test tubes

Total length (mm)	Usable length (mm)	Ø (mm)	Pack (u.)	Art. No.
200	85	25	1	233-0402-2
250	75	20	5	233-0402-8
250	75	15	5	233-0402-0
220	60	10	5	233-402-00

For pipettes

Total length (mm)	Usable length (mm)	Ø (mm)	Pack (u.)	Art. No.
500	85	7	5	233-0402-9
400	75	5	5	233-402-15

For further information, please contact scharlab@scharlab.com