



EMB plate 150mm

Precision balances EMB

KERN

Entry level laboratory balance with tremendous weighing performance.

- Ready for use: Batteries included, 9V, operating time approx. 12 h.
- AUTO-OFF function to preserve the battery, can be switched off
- Tare function facilitates formulation work
- Particularly flat design
- Secure and non-slip positioning with rubber feet

Technical data:

- Large LCD display, digit height 15mm
- Dimensions of weighing plate (plastic): 82mm, 105mm, 150mm
- Overall dimensions WxDxH 170x240x39mm
- Net weight approx. 0,5kg
- Permissible ambient temperature 5°C / 35°
- Power supply: 220-240 V

Available accessories like draft shield or mains adapter external are detailed below.

Optional DKD calibration certificate is available. It must be requested on the same order along with the balance.

Model	Readability (g)	Capacity (g)	Reproducibility (g)	Linearity (g)	Plate (mm)	Calibration	Optional cert.	Pack (u.)	Art. No.
EMB 100-3	0,001	100	0,001	±0,003	82	External	315-963127	1	00EMB100-3
EMB 200-2	0,01	200	0,01	±0,02	105	External	315-963127	1	00EMB200-2
EMB 200-3	0,001	200	0,001	±0,003	82	External	315-963127	1	00EMB200-3
EMB 600-2	0,01	600	0,01	±0,03	105	External	315-963127	1	00EMB600-2
EMB 1000-2	0,01	1.000	0,01	±0,03	150	External	315-963127	1	0EMB1000-2
EMB 2000-2	0,01	2.000	0,01	±0,03	150	External	315-963127	1	0EMB2000-2
EMB 500-1	0,1	500	0,1	±0,2	150	External	315-963127	1	00EMB500-1
EMB 1200-1	0,1	1.200	0,1	±0,3	150	External	315-963127	1	0EMB1200-1
EMB 2200-0	1	2.200	1	±2	150	External	315-963127	1	0EMB2200-0
EMB 3000-1	0,1	3.000	0,1	±0,3	150	External	315-963127	1	0EMB3000-1
EMB 6000-1	0,1	6.000	0,1	±0,3	150	External	315-963128	1	0EMB6000-1
EMB 5.2K1	1	5.200	1	±3	150	External	315-963128	1	315EMB52K1
EMB 5.2K5	5	5.200	5	±10	150	External	315-963128	1	315EMB52K5

Weighing elements
 in **page 790**



Visit our new web www.scharlab.com