


Proline® Plus BIOHIT adjustable-volume automatic pipettes
SARTORIUS

Sartorius Biohit's mechanical pipette family Proline Plus is designed to offer comfort and quality in every day manual pipetting. Proline Plus pipettes are an excellent choice for both experienced laboratory professionals and students. Ergonomic design with light pipetting forces, comfortable handle and finger support for reduced risk of strain injuries. Minimised risk of contamination with Safe-Cone Filters and full autoclavability. Highly durable with strengthened structure also for heavier use.

Channels	Volume range (µl)	Increment (µl)	Test volume (µl)	Inacc. (%)	Impr. (%)	Pack (u.)	Art. No.
1	0.1-3	0.002	3	1.30	0.80	1	0000728010
			1.5	2.40	1.60		
			0.3	10.00	6.00		
1	0.5-10	0.01	10	1.00	0.60	1	0000728020
			5	1.50	1.00		
			1	2.50	1.50		
1	2-20	0.02	20	0.90	0.40	1	0000728030
			10	1.20	1.00		
			2	3.00	2.00		
1	5-50	0.10	50	1.00	0.30	1	0000728040
			25	1.40	0.50		
			5	3.00	1.50		
1	10-100	0.10	100	0.80	0.15	1	0000728050
			50	1.00	0.40		
			10	2.00	1.00		
1	20-200	0.20	200	0.60	0.15	1	0000728060
			100	0.80	0.30		
			20	2.00	0.80		
1	100-1000	1.00	1000	0.60	0.20	1	0000728070
			500	0.60	0.20		
			100	1.00	0.40		
1	500-5000	10.0	5000	0.50	0.20	1	0000728080
			2500	0.60	0.30		
			500	2.00	0.60		
1	1-10 ml	20.0	10000	0.60	0.20	1	0000728090
			5000	1.20	0.30		
			1000	3.00	0.60		
8	0.5-10	0.01	10	1.50	1.00	1	0000728120
			5	2.50	2.50		
			1	4.00	4.00		
8	10-100	0.10	100	0.70	0.25	1	0000728130
			50	1.00	0.70		
			10	3.00	1.50		
8	30-300	0.20	300	0.60	0.25	1	0000728140
			150	1.00	0.50		
			30	2.00	1.00		
12	0.5-10	0.01	10	1.50	1.00	1	0000728220
			5	2.50	2.50		
			1	4.00	4.00		
12	10-100	0.10	100	0.70	0.25	1	0000728230
			50	1.00	0.70		
			10	3.00	1.50		
12	30-300	0.20	300	0.60	0.25	1	0000728240
			150	1.00	0.50		
			30	2.00	1.00		

For further information, please contact scharlab@scharlab.com