


eLINE® BIOHIT electronic pipettes
SARTORIUS

The electronic Biohit eLINE pipette family has been providing high precision and comfort for years. The comprehensive range of pipetting modes in eLINE reduces the need for several work stages and enables liquid dispensing twice as fast as a mechanical pipette. Proven accuracy and precision with DC-motor concept. Fully electronic operations guarantee user-independent results and prevent RSI. One eLINE covers the volume range of two mechanical pipettes - fewer pipettes needed. Several charging options. Liquid dispensing from the air, without touching the vessel wall, with eLINE 0.1-5µl model.

Channels	Volume range (µl)	Increment (µl)	Test volume (µl)	System. error (%)	Random error (%)	Pack (u.)	Art. No.
1	0.1-5	0.05	5	1.20	0.70	1	0000730011
			2.5	1.70	1.20		
			0.5	9.50	8.00		
1	0.2-10	0.05	10	0.90	0.40	1	0000730021
			5	1.00	0.70		
			1	3.50	2.30		
1	5-120	0.50	120	0.60	0.15	1	0000730041
			60	0.80	0.20		
			12	2.50	1.00		
1	10-300	1.00	300	0.40	0.15	1	0000730061
			150	0.60	0.20		
			30	1.50	0.80		
1	50-1000	5.00	1000	0.40	0.15	1	0000730081
			500	0.60	0.20		
			100	1.50	0.50		
1	100-5000	10.0	5000	0.50	0.15	1	0000730101
			2500	0.80	0.20		
			500	1.00	0.40		
8	0.2-10	0.05	10	0.90	0.50	1	0000730321
			5	1.50	0.80		
			1	4.00	3.00		
8	5-120	0.50	120	0.80	0.20	1	0000730341
			60	0.70	0.30		
			12	3.00	1.50		
8	10- 300	1.00	300	0.50	0.20	1	730361
			150	0.70	0.30		
			30	2.00	1.00		
8	50-1200	5.00	1200	0.50	0.20	1	0000730391
			600	1.00	0.30		
			120	2.50	1.00		
12	0.2-10	0.05	10	0.90	0.50	1	0000730421
			5	1.50	0.80		
			1	4.00	3.00		
12	5-120	0.50	120	0.80	0.20	1	0000730441
			60	0.70	0.30		
			12	3.00	1.50		
12	10- 300	1.00	300	0.50	0.20	1	0000730461
			150	0.70	0.30		
			30	2.00	1.00		
12	50-1200	5.00	1200	0.50	0.20	1	0000730491
			600	1.00	0.30		
			120	2.50	1.00		