

Cylinders, graduated

Graduated cylinders

All our graduated cylinders are calibrated “to contain”.

Graduated cylinders come with fused on, high contrast, ceramic enamel blue printed lines and inscriptions.

• Correct handling of graduated cylinders

Fill with the liquid until the meniscus is exactly at the desired ring mark. The bottom of the meniscus should be read at eye level. The inside cylinder wall should not be wetted above the mark.



Graduated cylinder, to contain, Class A, with pour spout, serialized and certified

Borosilicate glass. Marked with individual lot number. Supplied with certificate of performance. Hexagonal base. Marked with a blue single metric scale. Complies with ISO 4788.

Art. No.	Capacity (ml)	Division (ml)	Units
073-201750	10	1/5	2
073-201751	25	1/2	2
073-201752	50	1/2	2
073-201753	100	1/1	2
073-201754	250	2/1	2
073-001755	500	5/1	1
073-001756	1.000	10/1	1

Graduated cylinder, to contain, Class B, with pour spout

Soda glass. Hexagonal base. Marked with a blue single metric scale. Complies with ISO 4788.

Art. No.	Capacity (ml)	Division (ml)	Units
073-101743	10	1/10	2
073-101744	25	1/5	2
073-101745	50	1/2	2
073-101746	100	1/1	2
073-101747	250	2/1	2
073-101748	500	5/1	1
073-101749	1.000	10/1	1



Cylinders, graduated