

**BURETTES AND DISPENSERS: Digital burettes**

**Solarus® titration digital burettes**  
**HIRSCHMANN**

A titration unit with a reagent recirculation system for manual titration directly from the reagent bottle. Power supplied by a solar cell.

**Main features:**

- Media recirculating system for simple and rapid venting without reagent loss
- Easily replaceable valves
- Handwheels with soft touch surface
- Clearly legible display, regardless of the viewing angle
- A light protection window for the valve block is also included in the scope of delivery
- Individual calibration of user-specified adaptations with "Quick-Cal"
- Serial bidirectional interface (RS232 and USB)

Provided with 3 adaptors (GL32, GL38 and S40) suction hose, return tube, key, manual and certificate.

Description	Inaccuracy (R)	Reproducibility (CV)	Pack (u.)	Art. No.
Solarus® 10ml	0,3%	0,1%	1	6709392010
Solarus® 20ml	0,2%	0,1%	1	6709392020
Solarus® 50ml	0,2%	0,1%	1	6709392050

**Accessories. (Check for other accessories)**

Description	Pack (u.)	Art. No.
Liqui-Soft® software	1	6709564010
RS232 interface	1	6709564005
USB interface	1	6709564006
Discharge tube, spiral	1	6709313027

**HIRSCHMANN®**

# Simply the best

**pipetus® –**

the pipette filler

- With LC display, everything under control
- Intelligent rapid charging technology
- High-tech material for precision pipetting


**ceramus® –**

the bottle top dispenser

- Ceramic piston made of extremely hard high-purity aluminium oxide ceramic
- Precision instead of sense of proportion
- Perfect ergonomics


**solarus® –**

the first solar-powered digital burette

- Sensitive and precise
- The Hirschmann recirculation system
- Solar technology for sustainable performance


**rotarus® –**

the high-class peristaltic pump

- Automatic pump head and tube detection with RFID technology
- External TFT touchscreen with user-friendly menu
- Rapid tool-free changing of pump heads


**Hirschmann Laborgeräte GmbH & Co. KG**

Hauptstraße 7 - 15 · 74246 Eberstadt Germany · Fon +49 7134 511 0 · Fax +49 7134 511 990

www.hirschmannlab.com · info@hirschmannlab.com