

## MICROBIOLOGY: Sterile cutlery

### Sampling tools GOSSELIN SAS

Cutlery made of white resistant polystyrene designed for sampling of liquids and aseptic solids. Individually packaged in a sealed bag and sterilized by irradiation.



Description	Pack (u.)	Art. No.
Knife 200mm	1.000	630-0KNIFE
Fork 155mm	1.000	630-00FORK
Small spoon 155mm	1.000	630-000001
Spoon 165mm	1.000	630-0SPOON
Multicutlery 180mm, volume 4ml	500	120-982897

## MICROBIOLOGY: Sterile swabs

### Sterile swabs DELTALAB

Swab dimensions:  
Wood and plastic: 150x2,5mm (head ±5mm)  
Aluminium: 150x0,9mm (cabeza ±1,5mm)

Metallic ones are used for sampling of narrow places with difficult access due to its flexibility.



A: Wood, B: PS breakable, 3: Aluminium

Labeled tube

Individual bag

Bagged 2 by 2

Description	Presentation	Pack (u.)	Art. No.
Wooden body and cotton head	Individual bag	2000	027-300200
Wooden body and cotton head	Bagged 2 by 2	1000	027-300210
PS breakable body and cotton head	Individual bag	2000	027-300201
PS breakable body and cotton head	Bagged 2 by 2	1000	0030021101
Aluminium body and cotton head	Individual bag	2000	027-300203
Wooden body and cotton head	Labelled tube	500	027-300250
PS breakable body and cotton head	Labelled tube	500	027-300261
Aluminium body and cotton head	Labelled tube	500	027-300251

### Sterile swabs with media

Description	Presentation	Pack (u.)	Art. No.
Peptone water + Neutralizer (peptone 3ml + 10% neutralizer)	Labelled tube	100	114-000005
Ringer's + Neutralizer (2,5ml Ringer's, 10% neutralizer)	Labelled tube	25	114-000007
Amies medium. Breakable plastic body and synthetic head	Labelled tube	100	120-605629