



J.P. Selecta

Water and oil baths "Precisterm"

1



Adjustable temperature range from ambient +5°C up to 110°C, stability $\pm 1^\circ\text{C}$.

Adjustable temperature range from ambient +5°C up to 200°C, stability $\pm 2^\circ\text{C}$.

- Extensive range of sizes and accessories
- Stainless steel interior and external case

Features

Capacity; from 2 litres up to 45 litres.

Double skinned tank of aisi 310 stainless steel with an exterior case of aisi 304 stainless steel.

Heating elements embedded in incoloy heat and corrosion resistant stainless steel.

Control panel

Illuminated mains power switch.

Hydraulic thermostat temperature control, synchronised with the "on" indicator lamp.

Safety thermostat cut out indicator lamp.

Safety

Din 12877.2 Standard.

The unit will cut out if there is insufficient liquid or if the bath temperature exceeds the set temperature.

A manual reset adjustable temperature safety thermostat is placed at the back of the unit.

Models

Order No.	Maximum Temperature °C	Capacity litres	Height / Widht / Depth (usable) cm	Height / Widht / Depth (exterior) cm	Power W	weight kg
0196000137	110	2	15x15x13	20x20x27	390	3.5
0196000138	110	5	15x30x14	20x37x27	600	6
0196000140	110	12	15x30x30	20x40x42	1100	8
0196000141	110	20	15x48x30	20x58x42	1500	10
0196000142*	110	45	15x62x50	22x69x62	2400	18

*With built in drain tap.

